

ESSENTIALLY BETTER



Classiko

The robust care bed for safety and comfort



STIEGELMEYERGROUP





Classiko

The robust and comfortable care bed

Comfort and safety for the resident, alleviation for care staff and an individual design – the Classiko unites many strengths to give good value for money. This **robust care bed** has the advantage of durability and a range of adjustment options that make everyday life more agreeable.

The **four-section mattress base** adapts to suit many different needs. For example, by raising the back and thigh rests and tilting the bed into the optional reverse-Trendelenburg position, the resident is moved into a comfortable sitting position for eating or watching television. With a height

adjustment range from 40 to 80 cm, the Classiko makes getting in and out of bed easy while ensuring a **back-friendly working position** for care staff.


The motorised model of the bed can be operated effortlessly with a **selectively lockable handset**.

An electrically adjustable bed is not necessary or desirable in all care situations, however. For this reason, Stiegemeyer also offers a **mechanical version** of the Classiko. By dispensing with cables and manually removable elements, this model meets the

safety requirements in psychiatric care particularly well. With a high safe working load of 225 kg, it is also suitable for heavy residents.

Care staff will find the mechanical Classiko easy to operate. Height adjustments via the pedal with the help of a **hydraulic pump** can be made just as easily as the force-assisted adjustments of the back and thigh rests. The operator makes savings in energy costs and expenses for leakage current measurements. The Classiko offers visual appeal with a large selection of different **head and footboards** and **wood decors** to suit every taste.

External dimensions/height adjustment range:



ESSENTIALLY BETTER

High safe working load

- electrically adjustable: 175 kg
- mechanically adjustable: 225 kg





Optimum protection for the resident

Safety sides and practical equipment options

The safety sides on the Classiko can be adjusted to meet the individual protection needs of the resident. The standard version of the bed is equipped with **full-length safety sides** on both sides of the mattress base for optimum safety.

If the resident's mobility needs to be supported, then the **combinable split safety sides** are a good choice. When raised at the head end only, this leaves enough room for getting in and out of bed. The split safety side can also be combined with a full-length safety side on the other side of the bed.

Independent residents with good mobility will not need safety sides at all – in this case, the Classiko can be fitted with a **fixed side panel**, which gives the bed a particularly homelike appearance.

With its equipment and accessory options, the Classiko becomes a comfortable bed adapted to your requirements.

The four-section steel wire mesh mattress base is robust and flexible. In the motorised model, the **auto-contour setting** and the **sliding backrest** help protect the resident from shear forces.

For taller residents, we recommend **bed extensions** to 210 or 220 cm. Matching mattress extension pieces and extended safety side bars are also available.

The Classiko is easy to move and manoeuvre. Thanks to its **centrally locking castors**, care staff can lock or unlock all four castors at once by pressing just one brake pedal. The motorised model is operated using a **handset**, which shows all of the options as easily understandable pictograms. Individual functions can be blocked for the resident's safety, if required.



▲ The Classiko bed is fitted with a full-length safety side as standard.



▲ The 4-section wire mesh mattress base is both sturdy and flexible.

▼ The Classiko care bed is optionally available with combinable split safety sides.



▼ The bed extension provides extra comfort for tall residents.



▲ Bed adjustments on the mechanical model use a pedal and are hydraulically assisted.



▲ Selectively lockable handset for easier operation.

Decors

The Classiko can be adapted to suit any room with a variety of attractive decors to choose from. The selection ranges from the light Lindberg Oak to the

striking dark Bella Noce Choco. Our customer service would be happy to show you the Stiegelmeier colours from our latest Wodego collection.



Technical data / Features

Dimensions and weights	Mechanical	Electric
External dimensions (depending on the wooden surround)	approx. 208 x 101 cm	approx. 208 x 101 cm
Mattress base (mattress dimensions)	200 x 90 cm	200 x 90 cm
Backrest length compensation	approx. 10 cm	approx. 10 cm
Safe working load, electric	-	175 kg
Safe working load, mechanical	225 kg	-
Ground clearance	approx. 15 cm	approx. 15 cm
Total weight	approx. 130 kg	approx. 130 kg
Standard features		
4-section mattress base	✓	✓
Wire mesh mattress base	✓	✓
Ball-bearing castors Ø 125 mm	✓	✓
Permanently unlockable castors	✓	✓
4 central locking castors	✓	✓
Height adjustment range 40 to 80 cm	✓	✓
Backrest adjustment to 70°	✓	-
Flat, raised leg position to 11°	✓	✓
Auto contour setting	-	✓
DBfK backrest length compensation	-	✓
Handset with selective locking function	-	✓
Two-bar full-length safety sides	✓	✓
A wide range of wooden surrounds in various decor combinations	✓	✓
Argentum metal finish	✓	✓
Accessory equipment		
Tiltable to reverse-Trendelenburg position	0	0
Three-bar full-length safety sides	0	0
Integral castors Ø 125 mm	0	0
Double castors Ø 150 mm	0	0
Tiltable to Trendelenburg and/or reverse-Trendelenburg position	0	0
Bed extension 210 or 220 cm	0	0
Pillow shelf	0	0
Wall deflection rollers, horizontal	0	0
Wall deflection rollers, vertical, at head end	0	0
Fixed side panel, instead of integrated safety sides	0	0
Sleeve for attaching universal holder	0	0
Combinable safety sides, split on one or both sides	0	0

✓ = Standard feature
 0 = Available as an option
 - = Not available

Safety and comfort for psychiatric care

The mechanical model of the Classiko is perfectly suited to the high safety requirements in psychiatric care. Since it has **no cables** and its components are protected against removal, it prevents injuries to residents and damage to the bed. At the same time, the comfort and **homelike appearance** of the Classiko enhance the feeling of well-being.

For more information on use in **psychiatric care**, go to www.stiegelmeier.com



Custodians of values
 Shapers of the future

STIEGELMEYERGROUP

Custodians of values

Our company was founded in Herford in 1900 and is still a family-run enterprise today. Respect and a sense of responsibility towards our customers, employees and partners form the basis for our day-to-day work. Based on these values and our many years of experience, we offer our customers all over the world high-quality products, solutions and services to meet the needs of the care and healthcare sectors.

Shapers of the future

In order to cope with rapid changes in the care and healthcare sectors, we develop all our products ourselves. Creating innovations and visions of what care will be like in future are the goals we pursue. We maintain regular contact with our customers in order to support them in their daily tasks. Implementing tomorrow's trends today is a successful tradition of ours.

Our expertise

Hospital


Nursing home

Homecare

Lifestyle beds

Furniture systems

Service



The Stiegelmeier-Group is your international partner for hospitals and care facilities and, with the Burmeier brand, also for caring in the home. We support our customers with practice-oriented products, solutions and services that we are constantly developing further.

- Founded in 1900 and still family managed
- Over 1,200 employees worldwide
- All core products manufactured in proprietary factories
- Product development and production in Germany
- Service, training and process optimisation at customer's premises
- International certification and export to more than 60 countries
- Distribution companies located in Belgium, Finland, France, Netherlands, Poland and South Africa

Our brands

STIEGELMEYER

BURMEIER



STIEGELMEYERGROUP

Hauptsitz / Headquarters

Deutschland / Germany
Stieglmeyer GmbH & Co. KG
Ackerstraße 42
32051 Herford
Phone +49 (0) 5221 185 - 0
Fax +49 (0) 5221 185 - 252
info@stieglmeyer.com
www.stieglmeyer.com

France
Stieglmeyer SAS
8 Allée des Chevreuils
69380 Lissieu
Phone +33 (0) 4 7238 - 5336
Fax +33 (0) 4 7238 - 5481
contact@stieglmeyer.fr
www.stieglmeyer.fr

Südafrika / South Africa
Stieglmeyer Africa (PTY) Ltd
The Estuaries · Century City
Cape Town, 7441
Phone +27 (0) 21 555 - 1415
Fax +27 (0) 21 555 - 0639
info.southafrica@stieglmeyer.com
www.stieglmeyer.co.za

Nederland / Netherlands
Stieglmeyer BV
Peppelenbos 7
6662 WB Elst
Phone +31 (0) 345 581 - 881
info@stieglmeyer.nl
www.stieglmeyer.nl

Suomi / Finland
Stieglmeyer Oy
Ahjokatu 4 A
15800 Lahti
Phone +358 (0) 3 7536 - 320
Fax +358 (0) 3 7536 - 330
stieglmeyer@stieglmeyer.fi
www.stieglmeyer.fi

Polska / Poland
Stieglmeyer Sp. z o.o.
86-212 Stolno
Grubno 63
Phone +48 56 6771 - 400
Fax +48 56 6771 - 410
info.stolno@stieglmeyer.pl
www.stieglmeyer.pl



Finden Sie Ihren Ansprechpartner auf unserer Website:
Find your local partner on our website:
www.stieglmeyer.com/partner

stieglmeyer.com