

ESSENTIALLY BETTER



# Quado

The compact miracle

# Always in the right place

## The mobile bedside cabinet Quado

The new bedside cabinet Quado by Stieglmeyer is a real compact miracle. With its functions reduced to the most important ones it is particularly suitable for daily hospital routines. At the same time 2 open compartments on both sides and a drawer that can be opened to both sides offer enough storage space for the personal belongings of the patient. With its neutral colours Quado adapts well to room furnishings and every hospital bed. If the overbed table is not needed, it can be lowered easily and damped.



▲ With its neutral colours Quado adapts well to room furnishings and every hospital bed.



◀ Thanks to gas spring assistance, the Quado can effortlessly be adjusted to the desired height of 85 to 115 cm.



▲ The Quado's overbed table can be lowered easily and damped.




▲ The extendable bookrest under the overbed table top offers the patient high reading comfort, without having to clear the table (optional).

# Transport

## Bed with bedside cabinet



Quado can be attached to the head or footboard for a joint transport with Stieglmeyer beds.



The easy transport of the Quado ensures relief during the daily routine. It can be attached to the head or footboard of almost all current Stieglmeyer hospital beds without effort, and thus can be transported together with the bed when a patient is relocated. The joint transport saves care staff from having to carry and, if applicable, document the personal items. During a simple transport the Quado can be rolled along with the bed. During a transport across thresholds or obstacles the cabinet should be lifted using the bed's

height adjustment. Once the Quado has arrived at its destination, it can be safely parked next to the bed thanks to its castors with brakes.



The Quado bedside cabinets can be pushed together for storage to save space.


# Hygiene

## Reliable reprocessing

The Quado is machine washable and can also be easily cleaned manually. This ensures optimum hygiene in the hospital.



▲ The chassis ensures high stability and straightforward cleaning.



The straight-line overbed table-top is particularly easy to clean.



▲ The hygiene-friendly design of the Quado lends itself to easy and thorough reprocessing.

# Accessories

## Useful additions



▲ With a towel hook on the rail of the bedside cabinet the towel is always at hand for the patient.

▲ With the longer hook toiletry bags can also be stowed close to the patient.

▲ An attachable USB extension cable allows charging of mobile devices in reaching distance of the patient.



▶ The drink holder offers secure hold for glasses and bottles.



◀ The optional additional compartment is suitable for magazines and further personal items.



▶ ▲ The optional sleeve for telephones and multimedia systems allows action in accordance with the patient's needs and the rooms' features.



# Technical data / Equipment

Technical data	
<b>Dimensions and weights</b>	
External dimensions with overbed table (W x D x H)	92 x 44.6 x 89.6–119.6 cm
External dimensions, overbed table retracted (W x D; incl. chassis)	56.4 x 44.6 cm
External dimensions body (W x D)	28.6 x 44.6 cm
Overbed tabletop height	85–115 cm
Overbed tabletop dimensions	40.7 x 57.4 cm
Total weight	ca. 25 kg
<b>Load capacity</b>	
Overbed table (attached centrally)	max. 15 kg
Upper drawer pulled out (fixed centrally, tipping stability)	2,5 kg
Top surface	15 kg
<b>Equipment / Options</b>	
<b>Drawers</b>	
Small top drawer	✓
<b>Bookrest</b>	
	0

<b>Castors</b>	
65 mm double castor (4 castors with brakes)	✓
<b>Colours and materials</b>	
Body colour	Light grey NCS S1500-N
Body material	PP, PA
Chassis colour	Chrystal White NCS S0500N
Chassis material	Steel
<b>Miscellaneous</b>	
Lock for top drawer	0
Sleeves for telephone systems	0
Adaptor for extending the telephone sleeve	0
Suitable for automated cleaning and disinfection	✓
Bottle holder on top surface rail	0
Towel hook on rail	0
Additional compartment on the side	0
Long hook on rail	0
USB extension cable	0

✓ = standard feature  
0 = available as an option

## STIEGELMEYERGROUP



The Stieglmeyer Group is your international partner for hospitals and care facilities and, with the Burmeier brand, also for home care. We support our customers with practical products, solutions and services, which we constantly develop further.

- Founded in 1900 and still family-run today
- Employs 1,100 members of staff worldwide
- Production of all core products in our own factories
- Product development and manufacture in Germany
- Servicing, training and process optimisation at the customer's site
- International certification and export to more than 60 countries
- Sales companies in Belgium, Finland, France, the Netherlands, Poland and South Africa

### Our brands





## STIEGELMEYERGROUP

### Hauptsitz / Headquarters

---

*Deutschland / Germany*  
**Stieglmeyer GmbH & Co. KG**  
Ackerstraße 42  
32051 Herford  
Phone +49 (0) 5221 185 - 0  
Fax +49 (0) 5221 185 - 252  
info@stieglmeyer.com  
www.stieglmeyer.com

*France*  
**Stieglmeyer SAS**  
8 Allée des Chevreuils  
69380 Lissieu  
Phone +33 (0) 4 7238 - 5336  
Fax +33 (0) 4 7238 - 5481  
contact@stieglmeyer.fr  
www.stieglmeyer.fr

*South Africa*  
**Stieglmeyer Africa (PTY) Ltd**  
The Estuaries · Century City  
Cape Town, 7441  
Phone +27 (0) 21 555 - 1415  
Fax +27 (0) 21 555 - 0639  
info.southafrica@stieglmeyer.com  
www.stieglmeyer.co.za

*Nederland / Netherlands*  
**Stieglmeyer BV**  
Peppelenbos 7  
6662 WB Elst  
Phone +31 (0) 345 581 - 881  
info@stieglmeyer.nl  
www.stieglmeyer.nl

*Suomi / Finland*  
**Stieglmeyer Oy**  
Ahjokatu 4 A  
15800 Lahti  
Phone +358 (0) 3 7536 - 320  
Fax +358 (0) 3 7536 - 330  
stieglmeyer@stieglmeyer.fi  
www.stieglmeyer.fi

*Polska / Poland*  
**Stieglmeyer Sp. z o.o.**  
86-212 Stolno  
Grubno 63  
Phone +48 56 6771 - 400  
Fax +48 56 6771 - 410  
info.stolno@stieglmeyer.pl  
www.stieglmeyer.pl



Finden Sie Ihren Ansprechpartner auf unserer Website:  
*Find your local partner on our website:*  
[www.stieglmeyer.com/partner](http://www.stieglmeyer.com/partner)

[stieglmeyer.com](http://stieglmeyer.com)