



Ravello

Safe and comfortable transport chairs

Ravello

The chair range for safe transfer

When patients or residents cannot move about independently, getting to other rooms can often be a challenge. What is needed is a transport chair which offers a high level of safety and comfort and which does not require a great deal of effort from staff. A homelike appearance contributes to preserving the dignity of patients and residents. The transport chairs in the Ravello range perfectly satisfy all these requirements. Due to the large variety of models and practical accessories, they can be adapted to meet individual wishes and requirements.

ESSENTIALLY BETTER Model facts

- Armrests can be lowered to provide a comfortable sitting position
- Optional side rests provide additional support
- Headrest, fitted as standard, offers enhanced safety
- Concealed care items on the back of the chair
- Optionally equipped with easy to attach table (accessory)



Ravello curo

- The lower leg rest and the fixed footrest lift simultaneously to follow backrest adjustments



Ravello sano

- 2 directional castors
- Retractable footrest for easier access



Ravello sano-s

- 4 directional castors
- Separately attached footrest can be retracted for better access



Ravello plus-s

- Seat height adjustable between 51 and 80 cm using a lever mechanism



▲ The length-adjustable footrest moves with the mattress base when this is adjusted and provides firm support for the user's feet.

Ravello curo

Reliable support

If a patient or resident needs to be transported lying down, the backrest on the Ravello curo can easily be reclined to an angle of 43°. The Ravello curo has large, stable front disc wheels and directional castors with brakes at the rear. This makes the chair effortless to manoeuvre.



If the Ravello sano and sano-s are swivelled into the lying position, only the lower leg rest moves at the same time.

ESSENTIALLY BETTER
The "s" factor

The "s" on the Ravello sano-s and Ravello plus-s models stands for "simple movement". These models have 4 swivel castors with a diameter of 125 mm or 150 mm. These chairs are particularly easy to steer and can be manoeuvred in even the tightest of spaces without collisions. This protects patients and residents and keeps the chair and room furnishings safe from damage.



The height can be adjusted between 51 and 80 cm using a simple lever mechanism on both sides.

Ravello sano and sano-s

Easy to get in and out of

The Ravello sano differs from the curo model in that the footrest is mounted separately. The advantage of this design is the freer access it allows to the seat, which is particularly beneficial for people with walking difficulties. The footrest can be retracted if the patient or resident wishes to sit down or get up. This allows

people in need of care to walk right up to the chair's seat and to feel firm ground beneath their feet when they stand up. This model is available as the sano with 2 directional castors, and as the variant sano-s with 4 directional castors.

- ▼ The footrest can be fully retracted and extended using a simple mechanism.

The footrest is fixed underneath the chair as an independent element. ▶



Ravello plus-s

For demanding care

If residents or patients have difficulty staying on their feet for longer than a moment, care staff can also quickly reach their physical limits. The Ravello plus-s helps in many practical ways in this situation. Thanks to the seat's height adjustment and tilting

options, the chair helps users in any position and makes it considerably easier to get in and out. Users feel the seat surface as soon as they start to sit down and do not have the feeling they are losing their footing.

- ▼ The seat can be tilted forward to facilitate standing up and sitting down.

- ▼ Adjustable armrests and a folding footrest assure unhindered access.

- ▼ The 4 swivel castors are particularly large and make them easy to manoeuvre.





▲ The optional tray can simply be attached to the armrests.



▲ The lowerable armrests allow unobstructed access from the front and the side.



▲ If the user requires additional support whilst sitting, there is a special headrest with lateral support available.



▲ There is an optional mounting bracket for drainage bags.



▲ Infusion stands can be inserted on both sides using the sleeves supplied as standard.

Equipment and options

Useful additions

A substantial standard model and innovative accessories transform the Ravello into a transport chair which people are happy to spend time in. The Ravello also allows medical aids to be attached securely and discreetly. All Ravello models are available in a radiant cobalt blue. The upholstery is resistant to standard disinfectants used in hospitals (concentration according to the list of the German Society for Hygiene and Microbiology, DGHM).



▲ Optional side rests promote an upright sitting position.

Technical data

Ravello	curo	sano	sano-s	plus-s
Dimensions and weights				
Overall width	approx. 68 cm	approx. 68 cm	approx. 68 cm	approx. 73 cm
Overall height	approx. 126 cm	approx. 126 cm	approx. 126 cm	approx. 126 cm
Overall depth	approx. 103 cm	approx. 103 cm	approx. 103 cm	approx. 109 cm
Seat	approx. 52 x 52 cm	approx. 52 x 52 cm	approx. 52 x 52 cm	approx. 50 x 48 cm
Safe working load	200 kg (static), 135 kg (dynamic)			
Castors				
Directional castors Ø 125 mm, with brakes at rear	✓	✓	✓	-
Directional castors Ø 125 mm, with brakes at front	-	-	✓	-
Directional castors back and front Ø 150 mm with brakes	-	-	-	✓
Front disc wheels Ø 200 mm	✓	✓	-	-
Adjustment ranges				
Reclining backrest	43°	43°	43°	30°
Lowerable arm rests	✓	✓	✓	✓
Footrest, moves with the leg rest	✓	-	-	-
Footrest, retractable	-	✓	✓	-
Footrest, folds up	-	-	-	✓
Height adjustment and inclining seat part	-	-	-	✓
Additional features				
Infusion stand sleeves on both sides	✓	✓	✓	✓
Attachable table	o	o	o	o
Side rest	o	o	o	o
Headrest	o	o	o	o
Drainage bag holder	o	o	o	o
Colour				
Upholstery colour cobalt blue	✓	✓	✓	✓



Fabric



STIEGELMEYERGROUP



Solutions for people who care

The Stieglmeyer Group is Germany's leading manufacturer of medical beds for hospitals, nursing homes and homecare. With modern bedside cabinets, furniture and digital assistance systems, it is also an innovative provider of comprehensive solutions. The family-run company successfully develops offers for international markets in more than 100 countries.

Our strengths:

- Product development at our Herford headquarters in Germany – production in four modern factories in Germany and Poland
- Factory in China for the Chinese market
- Subsidiaries and sales partners in many countries
- Top-class service with reliable spare parts supply
- Training and process optimisation for our customers
- Sustainable action with strong social commitment

Our expertise

- Hospital
- Nursing home
- Homecare
- Lifestyle beds
- Furniture systems
- Service

Our values

- Experienced
- Human
- Reliable
- Solution-focused
- Open-minded
- Confident

1,000
EMPLOYEES
WORLDWIDE

11
LOCATIONS
in Germany, Poland,
the Netherlands, Belgium,
France, Finland,
South Africa and China

Tailor-made
OFFERS
for all
MARKETS

MORE THAN
100,000
beds
produced
every year

International
CERTIFI-
CATION

GERMAN
ENGINEERING

UP TO
15-YEAR
spare part availability

DISTRIBUTION
in more than
100
countries

GERMAN
ENGINEERING

Our brands





STIEGELMEYERGROUP

Solutions for people who care



Hauptsitz / Headquarters

Deutschland / Germany
Stieglmeyer GmbH & Co. KG
Ackerstraße 42
32051 Herford
Phone +49 (0) 5221 185-0
Fax +49 (0) 5221 185-100
info@stieglmeyer.com
www.stieglmeyer.com



Hauptsitz / Headquarters

Deutschland / Germany
Burmeier GmbH & Co. KG
Industriestraße 53
32120 Hiddenhausen
Phone +49 (0) 5223 9769 - 0
Fax +49 (0) 5223 9769 - 090
info@burmeier.com
www.burmeier.com

France
Stieglmeyer SAS
16 Avenue de l'Europe
Bâtiment SXB1
67300 Schiltigheim
Phone +33 (0) 4 7238 - 5336
contact@stieglmeyer.fr
www.stieglmeyer.fr

South Africa
Stieglmeyer Africa (PTY) Ltd.
The Estuaries - Century City
Cape Town, 7441
Phone +27 (0) 21 555 - 1415
Fax +27 (0) 21 555 - 0639
info.southafrica@stieglmeyer.com
www.stieglmeyer.co.za

Nederland / Netherlands
Stieglmeyer BV
Peppelenbos 7
6662 WB Elst
Phone +31 (0) 345 581 - 881
info@stieglmeyer.nl
www.stieglmeyer.nl

Suomi / Finland
Stieglmeyer Oy
Ahjokatu 4 A
15800 Lahti
Phone +358 (0) 3 7536 - 320
Fax +358 (0) 3 7536 - 330
stieglmeyer@stieglmeyer.fi
www.stieglmeyer.fi

Polska / Poland
Stieglmeyer Sp. z o.o.
86-212 Stolno
Grubno 63
Phone +48 56 6771 - 400
Fax +48 56 6771 - 410
info.stolno@stieglmeyer.pl
www.stieglmeyer.pl

China
Stieglmeyer China Co. Ltd.
Poly Center, No. 8 South Wenhua Road
Dongping New Town, Shunde District
528300 Foshan
Phone +86 757 29878710 - 0
info@stieglmeyer.cn
www.stieglmeyer.cn



stieglmeyer-group.com