

ESSENTIALLY BETTER



Accessories for hospital beds

Enhanced flexibility and safety

International Edition

STIEGELMEYERGROUP





Contents

- 4 Patient lifting poles
- 6 Infusion stands
- 8 Standard sliding rails
- 10 Universal holders
- 12 Secure hold for oxygen bottles
- 13 Precision fitting of extension devices
- 14 Accessories on the footboard
- 16 Bed extensions
- 18 Lamps
- 19 Handsets
- 20 Safety sides
- 21 Accessories overview
- 22 The Stiegmeyer-group

Accessories from Stiegmeyer for your customised dream bed

Better quality of life and less stress

Even the standard versions of Stiegmeyer's hospital beds offer a wide range of options. With the right accessories, they become tailor-made dream beds that cater to the patient's needs and considerably relieve the burden on care staff. The central focus is on optimal medical care, improving safety for the patient, and ensuring a higher quality of everyday life. The accessories can also be purchased retrospectively and easily attached or removed. Enhance your establishment's public image, improve the standard of patient care, and protect your valuable investments.



Patient lifting poles
for comfort and support

Page 4



Infusion stands
for a reliable supply

Page 6



Standard sliding rails
for intensive care

Page 8



Universal holders
and urine bottle holders

Page 10



Secure hold
for oxygen bottles

Page 12



Precision fitting
of extension devices

Page 13



Accessories on the footboard
Always up-to-date and at the ready

Page 14



Bed extensions
More room and comfort for patients

Page 16



Lamps
Directly on the bed

Page 18



Handsets
always within easy reach

Page 19



Safety sides
for maximum protection

Page 20



Accessories overview

Page 21

Patient lifting poles for comfort and support

In order to lie comfortably in bed for a lengthy period of time, everyone has to change their position every so often. However, this is often a challenge for physically restricted patients in the hospital – and puts a strain on the caregivers who help them. Our patient lifting pole supports patients in **getting themselves** into a comfortable position on their own.

The patient lifting pole is inserted into the adapter sleeves on the right or left side of the hospital bed at the head end. The oval

metal tube – made of stainless steel, of course, at Stieglmeyer – is equipped with an **ergonomically pleasing grab handle**. The patient can raise himself/herself with little effort. We also offer an optional patient lifting pole with a double hook for infusions. They allow safe, space-saving treatment of the patient.

Crutches should always be close at hand. Our practical **crutch holders** are simply attached to the head and footboards of the bed or to the railing of the bedside cabinet.



The patient lifting pole helps with comfortable positioning in the bed.



Double hooks for infusions ensure a reliable supply for the patient. The soft-touch grab handle is particularly non-slip.




▲ If extensions are to be fitted above the patient, then the traction frame is ideal for this purpose. It offers longitudinal and transverse bars which can be variably positioned above the bed



▲ The patient lifting pole can simply be inserted into the adapter sleeves at the head end on either side of the bed.



▲ Patient lifting pole on the Evario with a fixed headboard.



▲ The crutch holders can be attached to the bed or bedside cabinet.

Infusion stands for a reliable supply

Infusion stands must be flexible and offer great stability at the same time. Models attached directly to the bed have advantages – they cannot become a stumbling block for patients and visitors and can be positioned close to the patient.

Our height-adjustable infusion stands can be inserted into the adapter sleeves at the head or foot end and enable all round patient supply with their four hooks. The **stand with a 45° bend** prevents damage to the wall of the room when the bed is moved to the Trendelenburg position, for example. It carries a total load of 8 kg, a single hook can hold up to 3 kg. For the **vertical infusion stand**, hooks with a safe working load of 5 kg each are also available as an option. For even higher requirements – such as for rinsing – a **heavy-duty infusion pole** is available with 2 hooks, each carrying 10 kg.

If there is a patient lifting pole on the bed, a self-locking **infusion holder with 4 hooks** can be attached to it. For our **Mobilo stretcher**, a **double telescopic infusion stand** is available, which can be attached to all 4 corners and can be stowed away at the foot end when not in use using an optional holder. Thanks to its swivelling hooks, it takes up little space there. The **sturdy mounting bow** is particularly suitable for the intensive care unit, and medical equipment can be securely attached to it.



▲ The mounting bow is suitable for use in intensive care units.



▲ The vertical infusion stand is the classic standard solution.



▲ The infusion stand with a 45° bend prevents collisions when adjusting the bed.



▲ The heavy-duty infusion pole carries loads of up to 10 kg per hook.



▲ The infusion holder with 4 hooks can be attached to the patient lifting pole.



▲ The telescopic infusion stand for the Mobilo stretcher can be stowed away at the foot end in a space-saving way.



Stainless steel standard sliding rails are a practical solution for attaching intensive care equipment to universal beds.



The additional adaptable standard sliding rail provides sufficient clearance between the lowered mattress base and the floor.



◀ Sliding rail on the bed frame (head end).



◀ Sliding rail on the backrest.



◀ Sliding rail on the Protega safety side.

Stainless steel sliding rails

Stainless steel standard sliding rails are **extremely stable** and provide a secure hold for medical equipment and drainage systems on the hospital bed. They are suitable for use in all hospital units, especially in intensive care units.

The **hook-on stainless steel sliding rail** can be hung from the accessory holders on the side of the bed. The support on the **adaptable standard sliding rail** is bent upwards, so that the sliding rail maintains sufficient clearance from the floor even at the lowest position of the bed. The **stainless steel sliding rail for the infusion stand** is suitable for space-saving attachment of accessories at eye level – for example infusion pumps.

Our standard sliding rails enable, among other things, the positioning of drainage bags close to the patient. Practical drainage **tube holders** are available to secure the drainage tubes additionally. They ensure that the tubes cannot be bent and can be laid in an orderly way. Holders are available for attaching to the handle bars of the headboards or to the mobilisation support on the MultiFlex safety side.



▲ A tube holder is also available for the mobilisation support of the MultiFlex safety side.



▲ The standard sliding rail for the infusion stand enables accessories, such as infusion pumps, to be positioned at eye level in a space-saving way.



◀ The drainage tube holder can be attached to the handle bars of the head and footboards.



The universal accessory rail with plenty of room for accessories.



▲ The universal accessory holder on the Puro hospital bed with flexible attachment heights for the low-height adjustable mattress base



▲ The urine bottle holder for the Protega safety side on the Evario.



▲ The urine bottle holder for the Puro's MultiFlex[®] safety side system with flexible attachment heights for the low-height adjustable mattress base.



▲ The urine bottle holder for attaching to the side of the bed.

Universal and urine bottle holders

Accessories, drainage bottles and urine bottles need to be positioned safely along the sides of the hospital bed. Tubes must not be sharply bent or become detached, and liquids must not leak. In addition to our standard sliding rails, there are further individual solutions for providing a **secure attachment** and a **high standard of hygiene**.

The sturdy **universal accessory rail** offers plenty of room to hang several accessories next to each other. Depending on the bed model and safety side used, the **universal accessory holder** or the **universal accessory bow** can be attached flexibly to the side of the bed. They ensure that drainage bags do not lie unhygienically on the floor and that tubes are optimally protected.

Urine bottles should be easy to reach whenever necessary and reliably prevented from tipping over. Our **urine bottle holders** are ideal for meeting these requirements. We offer models with customised holders for a variety of hospital beds and equipment – for example, for MultiFlex[®] safety sides on our Puro, and Protega safety sides on the Evario.



The universal accessory bow holds drainage bags and tubes.



On the Puro, oxygen bottles of up to 5 litres can be attached to the holder on the stainless steel standard rail.

Secure hold for oxygen bottles

The reliable supply of oxygen to respiratory patients is always a top priority. Oxygen bottles must remain in the required position on the bed, even when the bed is moved. Our oxygen bottle holder enables bottles containing up to 2 litres to be flexibly positioned on the head and footboards. The holder is **simply hung** from the handle bar. It fits all head and footboard

versions, such as the Intercontinental head and footboard of the Evario or the Stelo head and footboard of the Puro.

The oxygen bottle holder that can be **hung from the stainless steel sliding rail** has an even higher safe working load: 5-litre bottles can be attached using this holder.

▼ On the Intercontinental head and footboard model, the oxygen bottle holder is hooked on underneath the handle bar.

▼ The stainless steel handle bar on the Rondo footboard provides a secure hold for the oxygen bottle holder.



The foot-end extension fitting on the Evario hospital bed provides ideal conditions for health-promoting therapy for the patient.

Precision fitting of extension devices

If extremities need to be positioned in trauma surgery units, using perforated rod systems, for example, one of the basic requirements is the ability to attach extension devices safely

and precisely to the hospital bed. To make this possible, an **exchangeable extension fitting** can be attached to the bed in place of the footboard – for example, to the Evario or Deka. Its two sturdy poles provide an ideal hold for the devices.



◀ The foot-end extension fitting can also be securely attached to the Deka hospital bed.



The writing surface holds the patient file in place to allow records to be filled in.



▲ The document holder in the modern Argentum colour goes well with every bed.

Accessories on the footboard always up-to-date and at the ready


Even in times of digitalisation, fast handwritten notes and entries are often the preferred alternative in everyday hospital life. **Our writing surface** provides the best conditions for this. It can be easily and flexibly attached to the head and footboards and folded down when not in use, thereby saving space. Paper is held securely in place during writing. The **large document holder** is used to store the documents. It also offers space for magazines and is hung on the handle bars of the head or footboards or the bedside cabinet.

Being able to address patients by name provides reassurance and an atmosphere of trust and is a help for doctors on their rounds. Our plastic name tag holder for a **name tag** is practical and hygienic.



▶ The plastic name tag holder can be positioned flexibly.

▼ A handy towel rail can be attached to the handle bars of the headboard and footboard.



Mattress base extensions

More room and comfort for patients

People in general are growing increasingly taller and consequently need more space in a hospital bed. To make best use of the integrated extension which our beds offer, we recommend the use of our tailor-made accessories.

The **mattress extension piece** fills the gap between the mattress and the footboard and ensures maximum comfort. Depending

on the model of bed concerned, it can be attached in two ways without tools – either to the bed frame using a **fixed holder** or to an extension of the lower leg rest. The latter ensures that the extended lower leg rest will move along with the rest of the bed when it is adjusted. This means that even tall patients can enjoy healthpromoting positions with raised lower legs (e. g. orthopaedic positioning).



▲ On the Seta pro and many other models, an optional extension piece is easy to attach to the lower leg rest.



▲ The optional lower leg rest extension, with its matching mattress extension piece, moves together with the mattress base whenever this is adjusted, allowing even tall patients to benefit from a raised leg position.



▲ On some beds, the mattress base is enlarged with an extension section, including a mattress extension piece, that does not move with the lower leg rest – shown here on the Deka model.



▲ Some beds, such as the Evario here, can be configured with an integrated extendible lower leg rest extension. A matching mattress extension piece is also available.

ESSENTIALLY BETTER

Mattress base extensions with strong advantages

The **integrated extension of the lower leg rest** is recommended if you need to use the bed quickly at any time. This is always ready for extending immediately and then only needs the mattress extension piece to be inserted (optionally available for Evario, Evario select, Puro, Sicuro tera).

Separate parts and mattress pieces are an affordable alternative for extending the mattress base. These do not need to be ordered for each bed right from the start but can be fitted later or exchanged among individual beds (available for Evario, Puro, Seta pro, Sicuro tera, Deka, Vida).

◀ The standard bed extension of the Evario can be extended effortlessly.

Lamps

directly on the bed

Reading lamps from Stieglmeyer provide a pleasant light – without disturbing anyone else in the room. The lamps can be attached to the infusion sleeves at the head end of the bed. An integrated adapter prevents inadvertent rotation of the lamp arm.

Both LED and halogen models of the **Sola reading lamp** are available. In both cases, the lamp head is made of shatterproof plastic. Vents and heat dissipating cylinders ensure that heat build-up remains low.



▲ The Sola lamp has an attractively shaped lamp head made of unbreakable plastic.

The **Stella LED reading lamp** combines high-quality design with a long service life. The lamp head is made of highly polished stainless steel and provides night-time orientation lighting. The Stella model is also available as a version for attaching to the bedside cabinet.



▲ Our Stella reading lamp provides pleasant, cosy lighting for patients – without disturbing other people in the room.

Handsets

always within easy reach



▲ There is also a tailor-made holder for the standard handset.

Handsets should always be within easy reach of patients and nursing staff. At the same time, they should be kept in a safe place on the bed to prevent any damage. Stieglmeyer hospital beds already offer practical fixing options as standard.

For even more flexibility, we recommend our **holder** for attaching to the handle bars on the head and footboards. There are tailor-made versions available for the standard handset and the innovative LCD handset.

The **handset holder with a flexible shaft** is also particularly convenient for the patient. On selected models, it can be attached to the side of the bed within easy reach and can also be positioned exactly in the desired position.

◀ With a handset holder on the headboard or footboard, the LCD handset can be placed in or out of immediate reach of the patient.



▲ With the flexible transmitter holder at the side of the bed, the handset is always within ideal reach of the patient.

Safety sides for maximum protection

Classic hospital beds without integrated safety sides are still in use in many cases. For this purpose, we offer our **full-length safety sides**, which can be easily attached and removed to ensure the necessary flexibility. The safety sides offer maximum protection and comply with current normative requirements. When installed on one side only, the safety side can also help with orientation at night.

Foam leather padding is available for many safety side systems for even more comfort and safety. This feels soft and pleasant to the touch and can be reprocessed hygienically. The padding can be used on both full-length and split safety sides to protect particularly restless patients.



▲ This infill piece is available to increase the protective height in the middle section of the Protega safety side.



The protective cover pleasantly cushions the four ► Protega safety side elements.

▼ Beds without any safety sides can be retrofitted with 3-bar, full-length safety sides at any time.



▼ The 4-bar safety side is suitable for use with particularly high anti-decubitus mattresses.



▲ The protective cover for the full-length safety side moves up and down with it when the bars are raised and lowered.



▲ The split MultiFlex+ safety side can also be fully encased in an upholstered cover.

Overview

The design of the accessories may vary depending on the bed model.

Page		Evario	Evario select	Puro	Seta pro	Seta pro junior	Deka	Vida	Sicuro tera	Mobilo stretcher
Patient lifting pole										
4	Patient lifting pole, round tube	✓	✓	✓	✓	✓	✓	-	✓	-
5	Patient lifting pole, oval tube, with infusion hooks	✓	-	✓	✓	-	✓	✓	✓	-
5	Traction frame	✓	✓	✓	✓	-	✓	-	✓	-
6	Mounting bar (for transducer)	✓	-	✓	✓	-	-	-	✓	-
5	Grab handle	✓	✓	✓	✓	-	✓	✓	✓	-
7	Infusion holder (plastic) for patient lifting pole	✓	✓	✓	✓	-	✓	✓	✓	-
5	Crutch holder	✓	✓	✓	✓	✓	✓	✓	✓	-
Infusion stand										
6	Infusion stand	✓	✓	✓	✓	✓	✓	✓	✓	-
7	Infusion stand with a 45° bend	✓	✓	✓	✓	-	✓	✓	✓	-
7	Mobilo infusion stand	-	-	-	-	-	-	-	-	✓
7	Heavy-duty infusion pole	✓	✓	✓	✓	-	✓	-	✓	-
Standard sliding rails										
9	Adaptable standard sliding rail	-	-	✓	-	-	-	-	-	-
8	Hook-on standard sliding rail	-	-	✓	✓	-	-	✓	-	-
9	Standard sliding rail for infusion stands	✓	✓	✓	✓	-	✓	✓	✓	-
8	Standard sliding rail on the bed frame (head end)	✓	-	✓	-	-	-	-	-	-
8	Sliding rail on the backrest	✓	-	✓	-	-	-	-	✓	-
8	Standard sliding rail on Protega safety side	✓	✓	-	-	-	-	-	✓	-
Universal holders										
9	Drainage tube holder (for routing tubes)	✓	✓	✓	✓	-	✓	-	✓	-
9	Tube holder on MultiFlex+ safety side	-	-	-	✓	-	-	✓	-	-
11	Universal accessory holder	✓	-	✓	✓	✓	✓	✓	-	-
11	Urine bottle holder for MultiFlex safety side	-	-	✓	✓	-	-	✓	-	-
10	Hook-on universal accessory rail	-	-	✓	✓	-	-	✓	-	-
11	Urine bottle holder with or without anti-twist safeguard	✓	-	✓	✓	✓	✓	✓	-	-
11	Urine bottle holder for MultiFlex safety side	-	-	✓	✓	-	-	✓	-	-
11	Urine bottle holder for Protega safety side	✓	✓	-	-	-	-	-	✓	-
Oxygen bottle holders										
12	Oxygen bottle holder	✓	✓	✓	✓	-	✓	-	✓	-
12	Oxygen bottle holder for standard sliding rail	-	-	✓	✓	-	-	-	✓	✓
Fittings for extension devices										
13	Foot-end extension fitting	✓	✓	✓	✓	-	✓	-	-	-

✓ = Standard - = Not available

It may not be possible to attach accessories to all headboards, footboards or safety side variants

Overview

The design of the accessories may vary depending on the bed model.

Page		Evario	Evario select	Puro	Seta pro	Seta pro junior	Deka	Vida	Sicuro tera	Mobilo stretcher
	Holders on the footboard									
15	Document holder	✓	✓	✓	✓	✓	✓	-	✓	-
14	Attachable writing surface	✓	✓	✓	✓	-	✓	-	✓	✓
15	Towel holder	✓	✓	✓	✓	✓	✓	-	✓	-
14	Name tag holder	✓	✓	✓	✓	✓	✓	-	✓	✓
	Bed extensions									
17	Bed extension fitting, including mattress extension piece	✓	-	✓	✓	-	✓	✓	✓	-
17	Mattress extension piece for integrated lower leg rest extension	✓	✓	✓	-	-	-	-	✓	-
	Lamps (for use in infusion stand sleeve)									
18	Sola LED reading lamp	✓	✓	✓	✓	-	✓	-	✓	-
18	Stella LED reading lamp	✓	✓	✓	✓	-	✓	-	✓	-
	Handset holders									
19	Handset holder for LCD handset	✓	-	✓	-	✓	-	-	-	-
19	Handset holder for standard handset	✓	-	✓	✓	✓	-	-	-	-
19	Handset holder with a flexible shaft	-	-	-	-	-	✓	✓	-	-
	Safety sides									
20	Safety side, can be retrofitted, 3 bars	-	-	-	-	-	✓	-	-	-
20	Safety side, retrofittable, 4 bars	-	-	-	-	-	✓	-	-	-
20	Foam leather cover for full-length safety sides	-	-	✓	✓	-	✓	-	-	-
20	Foam leather cover for MultiFlex safety side	-	-	✓	✓	-	-	✓	-	-
20	Foam leather cover for Protega safety side	✓	✓	-	-	-	-	-	✓	-
20	Centre infill piece for Protega safety side	✓	✓	-	-	-	-	-	✓	-

✓ = Standard - = Not available

It may not be possible to attach accessories to all headboards, footboards or safety side variants



Custodians of values Shapers of the future

STIEGELMEYERGROUP

Custodians of values

Our company was founded in Herford in 1900 and is still a family-run enterprise today. Respect and a sense of responsibility towards our customers, employees and partners form the basis for our day-to-day work. Based on these values and our many years of experience, we offer our customers all over the world high-quality products, solutions and services to meet the needs of the care and healthcare sectors.

Shapers of the future

In order to cope with rapid changes in the care and healthcare sectors, we develop all our products ourselves. Creating innovations and visions of what care will be like in future are the goals we pursue. We maintain regular contact with our customers in order to support them in their daily tasks. Implementing tomorrow's trends today is a successful tradition of ours.

Our expertise

Hospital

Nursing home

Homecare

Lifestyle beds

Contract furniture

Service



The Stiegelmeier-Group is your international partner for hospitals and care facilities and, with the Burmeier brand, also for home care. We support our customers with practical products, solutions and services, which we constantly develop further.

- Founded in 1900 and still family-run today
- Employs 1,100 members of staff worldwide
- Production of all core products in our own factories
- Product development and manufacture in Germany
- Services, training and process optimisation at the customer's site
- International certification and export to more than 60 countries
- Sales companies in Belgium, Finland, France, the Netherlands, Poland and South Africa

Our brands





STIEGELMEYERGROUP

Hauptsitz / Headquarters

Deutschland / Germany
Stieglmeyer GmbH & Co. KG
Ackerstraße 42
32051 Herford
Phone +49 (0) 5221 185 - 0
Fax +49 (0) 5221 185 - 252
info@stieglmeyer.com
www.stieglmeyer.com

France
Stieglmeyer SAS
16 Avenue de l'Europe · Bâtiment SXB1
67300 Schiltigheim
Phone +33 (0) 4 7238 - 5336
Fax +33 (0) 4 7238 - 5481
contact@stieglmeyer.fr
www.stieglmeyer.fr

South Africa
Stieglmeyer Africa (PTY) Ltd
The Estuaries · Century City
Cape Town, 7441
Phone +27 (0) 21 555 - 1415
Fax +27 (0) 21 555 - 0639
info.southafrica@stieglmeyer.com
www.stieglmeyer.co.za

Nederland / Netherlands
Stieglmeyer BV
Peppelenbos 7
6662 WB Elst
Phone +31 (0) 345 581 - 881
info@stieglmeyer.nl
www.stieglmeyer.nl

Suomi / Finland
Stieglmeyer Oy
Ahjokatu 4 A
15800 Lahti
Phone +358 (0) 3 7536 - 320
Fax +358 (0) 3 7536 - 330
stieglmeyer@stieglmeyer.fi
www.stieglmeyer.fi

Polska / Poland
Stieglmeyer Sp. z o.o.
86-212 Stolno
Grubno 63
Phone +48 56 6771 - 400
Fax +48 56 6771 - 410
info.stolno@stieglmeyer.pl
www.stieglmeyer.pl



Finden Sie Ihren Ansprechpartner auf unserer Website:
Find your local partner on our website:
www.stieglmeyer.com/partner

stieglmeyer.com