Gigant



Clinic | Care | **Homecare** | Assist





Gigant -

The care bed for heavier residents weighing up to 280 kg

Advantages

Quality and durability

The "Gigant" care bed is extremely sturdily constructed to stand up to the rigours of daily care. High-quality components make this a particularly durable, sturdy and low maintenance care bed that significantly reduces follow-on costs.

Safe working load of 350 kg

The mattress base of the Gigant care bed can accommodate a weight of up to 350 kg. Residents weighing up to 280 kg will appreciate the generous 120 cm wide mattress base, which can be extended to a length of 220 cm if necessary (bed extension available as an optional accessory).

Mobility and comfort

Early mobilisation of those in need of care is of increasing importance in the daily care routine of today. This is why the Gigant assists this process with special functions. The mobilisation support in the combinable split safety sides (optional) with its ergonomic knob provides the occupant of the bed with a secure grip and enhanced safety. The combinable split safety sides also comply with the requirements in place to restrict the use of restraint measures. It can also be combined with full-length safety sides.

Design details



Locating sleeve for mobilisation support (left) and locking box (right)



Emergency release of backrest



Extremely sturdy safety side construction



Plastic panels on the chassis for a more home-like feel

Gigant heavy-duty care bed

Features

- Mattress base 120 x 200 cm, one-piece, not dismantleable, mattress base height adjustable from 40 80 cm, backrest, thigh rest and mattress base height are adjustable with electric motors
- Reverse-Trendelenburg position can be set by tilting the mattress base (using the handset)
- Emergency release of backrest
- · Metal parts finished with polyester powder coating, colour argentum
- Wooden surround in beech decor/solid beech
- Safe working load: 350 kg
- Max. patient weight: 280 kg

Manufactured to: standard IEC 60601-2-52:2010



Accessories

(available at extra cost)



Clamp-on extra-high safety side rail



Combinable split safety sides



Reading lamps



Bed extension including longer safety sides and side boards



Foam leather cover

Order numbers and external dimensions

59.6631.1 **Gigant** heavy-duty care bed, length 222 cm, width 135 cm

19.1400.1 Patient lifting pole with triangular handle (included in the supply package)

Burmeier GmbH & Co. KG

Pivitsheider Straße 270, 32791 Lage, Germany

Phone: +49 (0) 5232 9841 - 0 Fax: +49 (0) 5232 9841 - 41 Email: info@burmeier.com www.burmeier.com

