





Lippe IV

Dimensional sketch


Bed dimensions [mm]

			
Model	80cm	90cm	100cm
Max. width B1	785	885	985
Width of mattress frame B2	750	850	950
Chassis width B3	690-1100		
Model	190cm		200cm
Bed length L	1850		1950
Chassis length L1	1195		
Distance Chassis to head end L2	369		369
Distance Chassis to foot end L3	286		286
Width for ground clearance L4	930		
Ground clearance H	100 (optional +70)		
Bed extension	Head and foot side each 100		
Adjustment range of the lying surface height	290 - 740		
Tilt in low foot position	-		



Lippe IV

Dimensional sketch



We reserve the right to make technical changes. Slight colour deviations of the originals are due to printing technology. Status 05/2024.