

ESSENTIALLY BETTER



Beds and furniture for rehabilitation centres



STIEGELMEYERGROUP



STIEGELMEYER



The right bed
for every phase

Optimal care in the six phases of rehabilitation

In science, neurological rehabilitation is often described as a sequence of six phases. Stieglmeyer accompanies patients and carers through the entire process with tailor-made support:



A gentle return to independent living

Beds and furniture for all rehabilitation requirements

During successful rehabilitation, a patient goes through several phases with different needs. These range from a hospital bed in the intensive care unit to a versatile care bed and then a lifestyle bed with supporting functions that can also be used at home.

For these phases, Stieglmeyer offers all the right products from a single source. Our extensive experience in all areas of care results in beds and furniture that meet individual requirements and complement each other perfectly.

The recovery and mobilisation of the patient is aided by split safety sides, effortless height adjustment for easy entry and exit, and flexible mattress bases and mattresses for health-promoting positioning. The homelike design of our beds, bedside cabinets and wardrobes is also beneficial since it helps people to feel at home at each stage of their rehabilitation and gives them new vitality.

Four divisions, one goal: our product range for rehabilitation

- Hospital
- Nursing Home
- Lifestyle Beds
- Contract Furniture

Contents

4	Libra	13	Carino	17	Conturo
6	Elvido vervo	14	Lindeo	18	Mattresses
8	Rehabilitation beds with Vario Safe	15	Quado	20	Contract furniture
10	Quilted headboards and footboards	16	Vitano	24	Tables and chairs
12	Puro	17	Lumano	26	The Stieglmeyer-Group



Vario Safe Libra

Specially designed to meet a wide range of individual needs

Patients in rehabilitation need the best possible protection against falls. At the same time, however, **unrestricted mobility** remains decisive when it comes to restoring their physical constitution. The Libra bed, a **low-height bed** with optional **split safety sides**, offers both a high level of safety and full freedom of movement. Its lowest position is close to the floor, at around just 25 cm, which minimises the risk of fall injuries, while the raised split safety sides at the head end still leave enough space for getting in and out of bed. The model is also available as a **Libra partner double bed**, so that spouses or partners can stay overnight with the patient and contribute to the success of the rehabilitation.

- Large **height adjustment range**, approx. 25 – 80 cm, for fall prevention and back-friendly care
- High **safe working load** of 225 kg
- Optional **split safety sides** for individual protection in any situation
- **Homelike appearance** thanks to a choice of headboards and footboards, wood decors and upholstery
- **Vario Safe system** for tool-free replacement of safety side elements, headboards, footboards and panels
- Easy-to-use **bed extension** available for single and double beds



- ▲ Fixed upholstery with attractive fabrics and faux leathers underline the homelike impression.
- ◀◀ The Libra partner double bed offers all the care benefits of the single bed while allowing the patient's partner to stay overnight in the rehabilitation facility.
- ◀ With the flexible Vario Safe system it is possible to create tailored combinations of safety side elements, panels, headboards and footboards.



ESSENTIALLY BETTER

The homelike bed for a new start

The Libra is ideal for the further, medical and long-term rehabilitation phases (C, D and F). The optional Vario Safe system makes it possible to use the safety sides precisely as required to suit the patient's abilities and needs. The elegant, homelike appearance of the bed reinforces the feeling of returning to a normal and self-determined life.

Vario Safe Elvido vervo

The versatile low-height bed for the discerning

The Elvido vervo has a large **height adjustment range** and is therefore also suitable for intensive care situations. It offers the best possible protection from fall injuries and an optimal working height, even for tall care staff. Thanks to numerous optional features, this model can become a personalised, tailor-made bed.

- Height adjustment range: approx. 25 – 80 cm
- All residents, no matter how tall or short, have secure floor contact when leaving the bed; **orientation and balance** are heightened thanks to tactile perception by the soles of the feet
- **OpenBus drive system**, a standard feature, allows the use of the optional Out-of-Bed system
- Option of split safety side with **Vario Safe system**
- Integrated **bed extension**, a standard feature, offers plenty of room for tall patients



- ▲ Time for a thorough foot care session? With the Vario Safe system, the footboard can be removed in seconds.
- ◀◀ With the fabric-upholstered head and footboards of the Nobla range, the Elvido vervo looks like an elegant designer bed.
- ◀ With many wood finishes to choose from, the Elvido vervo matches any room decor.

Phase C Phase D Phase F

ESSENTIALLY BETTER

Protection for patients at risk of falling


In the further, medical and long-term rehabilitation phases (C, D and F), the Elvido vervo is particularly effective in supporting patients who need a high degree of protection. Its optional Out-of-Bed system notifies carers if a patient at risk of falling does not return to bed. Its lowest position, at around 25 cm, and the optional split safety sides also protect patients without restricting their mobility or hindering their rehabilitation.



Rehabilitation beds with Vario Safe – the best solution for standard-compliant safety

If you choose the **Elvido vervo** or **Libra** models with split safety sides or without safety sides, you will enjoy using our flexible Vario Safe system. The Vario Safe principle is as impressive as it is easy: The safety side elements, the headboard and footboard as well as the side panels can be attached and removed **without tools**. This is possible thanks to **easy-to-use locks** that snap into place securely and cannot be accidentally opened by the patient.

The headboard, footboard and panels can be supplied on request with attractive **fabric or faux leather upholstery**. Beds with Vario Safe can be adapted to the **changing needs** of patients in just a few easy steps. Rehabilitation facilities can order their beds with basic equipment and enhance them by adding more parts at any time.



Replace or remove footboards as you wish

- A removable footboard allows better access to the patient, e.g. for foot care.
- Choice of high footboards for an increased feeling of well-being, or low footboards for a free view
- Choice of attractive wood decors, fabrics and faux leathers

! Safety sides can be freely combined: 110/90 or 90/110

Modular design and adjustable without tools

Replace headboards in seconds whenever required

- Homelike design variants, with or without handle bars
- Choice of attractive wood decors, fabrics and faux leathers
- Headboards, footboards and side panels can be in different colours

ESSENTIALLY BETTER

Vario Safe For added flexibility and safety

Beds with Vario Safe meet the changing protection needs of the patients, and their look can be customised to suit any taste. The flexible safety side provides standard-compliant protection without restricting mobility during rehabilitation.



◀ All Vario Safe elements are inserted into the bed frame and snapped into place with easy-to-use locking devices.



◀ To release the locks again, pull the orange safety catch towards you and press down on the lock.


Attach or remove split safety sides as required

- Two elements with lengths of 90 and 110 cm can be freely combined or used separately
- Use of only 110 cm at the head end provides standard-compliant safety side protection
- Use of only 90 cm safety side element at the head end gives a feeling of security and orientation and leaves plenty of room for getting in and out of bed

Headboards and footboards with Vario Safe




Decors with Vario Safe



Information about other fabrics and faux leathers is available on request.

Two quilted patterns are available for homelike upholstery





Homelike design

Elegant upholstery with optional quilting for our Nobla head/footboard range

A pleasant atmosphere like that of a good hotel room or a private bedroom enhances the positive effect of rehabilitation. The optional upholstered surrounds of our Elvido vervo and Libra care beds have a part to play here. And with an elegant quilting pattern, they look even cosier.

There are two patterns to choose from: The tile-shaped quilting, with horizontal and vertical lines, has a modern look and radiates


fresh energy. The classic diamond-shaped quilting is both stylish and homelike. The wide ranging colours in Stieglmeyer's fabric collection bring out the best of both patterns and can be combined to suit every taste. With the help of the Vario Safe system, the headboards, footboards and panels can be replaced whenever a new look is desired.

Two quilted patterns are available:



Diamond-shaped quilting

Diamond patterns on furniture and clothing are often associated with a British lifestyle. They look both classy and cosy.



Tile-shaped quilting

The rectangular quilting of the tile design looks lively and energetic – the best prerequisites for rehabilitation.



Headboard: Nobla, footboard: Nobla reha

▲ The classic diamond-shaped quilting creates a distinguished atmosphere.



Headboard: Nobla, footboard: Nobla reha

▲ The tile-shaped quilting has a modern, clear-cut look.

Puro

The bed for fast track recovery

The Puro combines all the advantages that make life easier for patients and nursing staff: a large height adjustment range for conveniently getting in and out of bed and back-friendly working, a modern **LCD handset** and the choice of a full-length or split **MultiFlex⁺** safety side, which is easy to operate in even the most hectic situations.

- **A high safe working load** makes the bed equally suitable for heavier patients weighing up to 250 kg
- **Integrated 28 cm bed extension** for tall people
- **Easy cleaning and maintenance** due to large-surface design and good access to all components

Phase
A

Phase
B

ESSENTIALLY BETTER

Optimal start to rehabilitation

The modern Puro hospital bed meets the exacting demands of acute treatment and early rehabilitation (phases A and B). Its LCD handset with only 3 buttons and a large backlit display can be operated even with limited motor skills. From the first moment of treatment onwards, the MultiFlex⁺ split safety side combines protection with the possibility of early mobilisation.

► In the comfortable sitting position, patients can once again participate in everyday life at eye level. Breathing is made easier and at mealtimes the risk of aspiration is reduced.

►► The flat surfaces of the Puro promote safe hygiene in both manual and automatic reprocessing.



◀ The optional LCD handset offers customised control levels for patients, care staff and technicians respectively. With the EasyCare function, care staff have quick access to all of the bed functions.

EasyCare
gives care staff
access to
all options



The bed for
parent-child
health cures

▲ An inviting combination: the Carino with custom-made furniture from Stieglmeyer.

Carino

The homelike bed with fresh ideas

Stable, durable, comfortable and extremely versatile – the Carino is a bed with **many possible applications**. Whether in rehabilitation facilities, psychosomatic treatment centres or residential homes, it ensures excellent quality of life. The Carino is particularly appreciated for its **very good price-performance ratio** and individual design possibilities.

With its robust metal mattress base, the bed easily withstands numerous cycles of use. Mattress base **widths of 90 and 100 cm** are available and a bed extension is available for taller users. The Carino's **manually adjustable backrest** ensures a comfortable position when reading or watching TV. The optional reading lamp provides pleasant lighting – fixing sleeves are already in place and can also be used for fixing a patient lifting pole. Another practical accessory is the space-saving side tray, providing convenient space for a book or mobile phone.

The Carino has a sleek and elegant design and can be adapted to your own personal taste with a wide range of available **headboard and footboard forms** and wood decors to choose from.

The textile **Softcovers** for the headboard, footboard and sides are an especially cosy design feature of the bed. Fabrics in attractive colours and patterns can be used to transform the Carino into the centrepiece of the room. With an optional fixed board along the side of the bed, the Carino can also function as a comfortable sofa during the day.

Phase
E

ESSENTIALLY BETTER

The bed with many application options

Our Carino serves many purposes in medical and care facilities. It promotes follow-up rehabilitation (phase E) thanks to its adjustable backrest, serves as a comfortable bed for health resort guests and is also ideal for staff living areas.



For use at
home too

The four-section mattress base allows adjustment to a variety of comfortable and health-promoting positions.



▲ With its space-saving chassis, the Lindeo can be placed into practically any existing bed frame.



▲ Making the bed is much easier with the mattress base at its highest position.

Lindeo

The comfort frame for almost any bed

The Lindeo deluxe frame is simply inserted into an **existing bed** frame and offers all the advantages of an electrically adjustable rehabilitation bed. This makes it the ideal choice for use at home or in facilities that attach particular importance to a warm and cosy ambience.

With its **lift function** from 29 to 74 cm, it offers everyone an ideal height for easily getting in and out of bed. At the maximum height, changing the sheets is easy and the mattress can be used as practical surface for packing a suitcase – ideal for a self-regulated return to everyday life.

The **four-section mattress base** can be adjusted to many comfortable and health-enhancing positions with the help of a handset. With the backrest and thigh rest raised, it creates an ideal position for relaxed reading and watching TV or for breakfast in bed.

Phase
E

ESSENTIALLY BETTER

Support in all phases of life

The Lindeo deluxe frame offers the advantages of long-lasting, gentle support in follow-up rehabilitation (phase E) and in all subsequent stages of life. It promotes mobility and ensures healthy sleep.



The space miracle



▲ With a built-in refrigerator and lockable drawer (optional features), the bedside cabinet offers a special level of comfort for the patient.

Quado

Mobile and space-saving

The Quado is a space-saving server with many comfort features. Its **drawer** and the **two roomy compartments** can be accessed from both sides and offer generous storage space for personal belongings. With its **height-adjustable body** and the **fold-down overbed table** in light grey, the server is ideal for meals taken in bed. Its good manoeuvrability and light-weight construction also make the Quado suitable for use as a trolley for medicines and care utensils. Thanks to its **large, smooth surfaces**, it has excellent hygiene properties.

Vitano

The hospital bedside cabinet

The Vitano is a **spacious**, easy-to-use and attractive bedside cabinet designed for use in hospitals. Its combination of an (optionally lockable) drawer, an open compartment and a cupboard provides space for many personal items. These can be stored where they are **easily accessible** and protected from unauthorised access. Thanks to the **rotatable cabinet body**, the front of the Vitano is always within easy reach of the patient. With its height-adjustable body and fold-down overbed table, the server is ideal for meals taken in bed.



▲ Several Quado units can be pushed into one another to save space when in storage.



▲ The Quado can be effortlessly adjusted from a height of 85 to 115 cm, thanks to its gas spring support.



▲ The cup holder (optional) provides good support for glasses and bottles.



▲ Modern connectivity, with optional holders for phones and multimedia systems.



With slide-out table

▲ The Lumano is available with various decors and handles.

Lumano

The practical bedside cabinet

The Lumano, with its particularly large top surface, increases the level of comfort in daily rehabilitation routines. It has an **integrated slide out table** under the top surface which almost doubles the available space.

The interior of the Lumano also leaves nothing to be desired: The drawers glide on **roller slides** and are therefore particularly **smooth-running and quiet**. A simple release mechanism allows the drawers to be easily removed for thorough cleaning.



▲ The Lumano can also be fitted with a fold-down overbed table.



Accessible from both sides

▲ The Conturo's robust overbed table can be easily adjusted by hand.

Conturo

The robust bedside cabinet

Thanks to its **robust design**, the Conturo is strong enough to withstand all the demands of daily rehabilitation routines. It also offers many additional benefits: All drawers and the storage compartment can be used **from both sides**.

The manual **height adjustment mechanism of the overbed table** is easy to operate. If the patient wants to use the overbed table as a book rest, the table can be tilted.



▲ The colour of the mechanism matches the decor of the Conturo.



In the Fowler position, Sky-therm® ensures significantly more even pressure distribution than cold foam.



In the Fowler position, the soft-therm mattress distributes pressure considerably more evenly to the back, buttocks and heels.

For more information, see the mattresses brochure.

Z1 – Sky-therm® mattress with wavy cut

For relaxed and healthy sleep, it is crucial to distribute pressure **evenly** over the body and to reduce pressure peaks. Many types of viscoelastic foam achieve good results here, but at the same time their low flexibility reduces the mobility of the patient. The user sinks into the layer of foam and is no longer able to move effortlessly. Our Sky-therm® mattress offers a convincing solution here – and at the same time adapts very well to every adjustment of the mattress base.

Sky-therm® is a foam that ensures even pressure distribution and also facilitates users' movement thanks to its **excellent resilience**. With Sky-therm®, the skin remains well supplied with blood and is in full contact with the flexible mattress surface. A permanent exchange of heat and moisture can then take place.

The **mobility** and well-being of the user are promoted when it is easier to move in bed and also to get up. In hospitals, early patient mobilisation has a major impact on recovery.



▲ The incisions in the Sky-therm® layer improve the point elasticity and ventilation of the mattress.



- Reinforced edge zones
- Excellent point elasticity due to the lattice of the Sky-therm® layer
- The wavy cut protects problem areas of the body and increases mattress flexibility

Z2 – Soft-therm mattress with fine perforations

The soft-therm material is similar to conventional visco foam in its pressure-distributing effect and also has the benefit of very high **elasticity** and tear resistance. Our soft-therm mattress contributes to effective pressure ulcer prevention without restricting the user's freedom of movement. At the same time, the material is breathable and creates a good microclimate.

The soft-therm mattress follows every adjustment of the mattress base with its particularly **high flexibility**. This is due not just to the elastic top layer but also to the cold foam base, which is designed with perforations to allow greater movement.



- Highly elastic soft-therm layer
- Fine perforations in the base layer enhance the flexibility of the mattress and improve the microclimate



ESSENTIALLY BETTER Optimal positioning for every care situation

During rehabilitation, the many adjustment options of our beds are often used to their full extent. The precise adjustment of the mattress base makes it possible to find an ideal health-promoting position for each patient. It is very important, therefore, that the adjustments are not thwarted by a rigid mattress. Our Z1 and Z2 mattresses adapt to the contour of the mattress base with high precision in order to promote recovery.

Contract furniture






For quality of life

Our contract furniture transforms every room into a home. Wardrobes, cupboards and desks come together to create a practical and elegant ensemble. Our many **attractive wood decors** and our **range of designs** allow ample scope for creating a personal style. With its rounded corner profiles, ergonomic handles and discreet lighting, this sturdy furniture is ideal for use in rehabilitation facilities. Patients can put their clothes and valuables in spacious wardrobes and drawers which keep them safe but still within easy reach.

ESSENTIALLY BETTER

Everything from one source

Based on a comprehensive analysis, we design and build completely individual service rooms, reception areas etc. for you:

-  Arrange an appointment
-  Room analysis on site
-  Evaluation and design
-  Customised construction
-  Precise on-site installation

Tailor-made and professional – and all from a single source.

Special designs possible

To be able to work as effectively as possible in any room, the available space has to be utilised to its fullest extent. We produce furniture that is precisely tailored to your rooms and to suit your requirements and wishes.



Wall-hung furniture provides better access for wheelchair users and makes it easier to clean the floor.



The practical, beautifully designed hallway furniture immediately creates a welcoming feeling.

Wardrobes with individually configurable compartments and drawers provide easily accessible storage space.



Refrigerators with ventilation systems can be built into various items of furniture.



▲ Individual cupboard divisions can be selected for all chests of drawers and wardrobes.



▲ Elegance and flair for private-option patient rooms, in the style of a luxury hotel.



◀ Beds and furniture in psychiatric facilities have to meet the highest safety standards. Sturdy constructions, well thought-out details and the absence of cables prevent patients from endangering themselves and others.



▲ Mobile wardrobe interiors for a quick and easy change of room or ward.



▲ Storage space and a homelike feel in the rehabilitation room: this cupboard unit with Havana Cherry and Chalk decors includes an integrated refrigerator, a suitcase stand and indirect lighting elements.



▲ Duty rooms as individual as the requirements of residents and staff. We perform a room analysis and use this to develop a concept that is specially tailored to you.



▲ A large drawer with pull-out desk top becomes a sturdy work surface in the service room.

Customised contract furniture for any type of room

Whether you have a duty room, residents' room or dining room in mind – with beds and contract furniture from Stieglmeyer, many interiors in your residential home can be fully furnished. We also manufacture shelving, medicine cabinets and reception desks with millimetre precision to suit your needs. Together with you, our specialists develop ideal room concepts in which all those involved are happy to live and work.

ESSENTIALLY BETTER Suitable for every room

With furniture from Stieglmeyer, you have many options. You can put together suitable room furnishings from a large selection of preconfigured elements. We can also make customised furniture for you, accurate to the millimetre, on request. Even small or awkwardly laid out rooms unfold their full potential with our contract furniture.



◀ Open and variable design for day care facilities: comfortable seating, spacious wardrobes and secure lockers.



Tables and chairs with character and flair

Whether in the reception area, the dining room, the patient's living area or the duty room: the seating is often the calling card of the entire facility. Stiegmeyer offers comfortable and attractive chairs and sofas in **various design lines**. **Ergonomic dimensions**, clear lines and elegant upholstery create an inviting ambience for a relaxed exchange at eye level. Discover our modern Mavo and Fena chair ranges and the exclusive armchairs and sofas of the Luna series.

The chairs are complemented by **matching tables**. Four-legged dining tables, column tables with bistro flair, and low side tables are available to meet every requirement.

ESSENTIALLY BETTER

Practical and versatile

Our Mavo and Fena chair series are robust and easy to clean, and are suitable for many purposes. As stacking chairs, they can be quickly moved aside if the room needs to be cleaned or used for an event. The heavy-duty models support a safe working load of up to 180 kg. Optional handle bars along the backrests make the chairs easier to carry.

The Mavo chair series is characterised by distinctive contours. Handle bars on the backrests make the chairs easier to stack. ▶



Playful curves: the Fena chair series and the column table create a cosy bistro atmosphere that will awaken a desire to return to everyday life.



◀ Comfortable sofas and wing chairs – the Luna seating range exudes a cosy authenticity.

Custodians of values Shapers of the future

Custodians of values

Our company was founded in Herford in 1900 and is still a family-run enterprise today. Respect and a sense of responsibility towards our customers, employees and partners form the basis for our day-to-day work. Based on these values and our many years of experience, we offer our customers all over the world high-quality products, solutions and services to meet the needs of the care and healthcare sectors.

Shapers of the future

In order to cope with rapid changes in the care and healthcare sectors, we develop all our products ourselves. Creating innovations and visions of what care will be like in future are the goals we pursue. We maintain regular contact with our customers in order to support them in their daily tasks. Implementing tomorrow's trends today is a successful tradition of ours.

Our expertise

Hospital

Nursing home

Homecare

Lifestyle beds

Contract furniture

Service



The Stiegmeyer-Group is your international partner for hospitals and care facilities and, with the Burmeier brand, also for caring in the home. We support our customers with practice-oriented products, solutions and services that we are constantly developing further.

- Founded in 1900 and still family managed
- Over 1,100 employees worldwide
- All core products manufactured in proprietary factories
- Product development and production in Germany
- Service, training and process optimisation at customer's premises
- International certification and export to more than 60 countries
- Distribution companies located in Belgium, Finland, France, Netherlands, Poland and South Africa

Our brands





STIEGELMEYERGROUP

Hauptsitz / Headquarters

Deutschland / Germany
Stieglmeyer GmbH & Co. KG
Ackerstraße 42
32051 Herford
Phone +49 (0) 5221 185 - 0
Fax +49 (0) 5221 185 - 252
info@stieglmeyer.com
www.stieglmeyer.com

France
Stieglmeyer SAS
16 Avenue de l'Europe
Bâtiment SXB1
67300 Schiltigheim
Phone +33 (0) 4 7238 - 5336
contact@stieglmeyer.fr
www.stieglmeyer.fr

South Africa
Stieglmeyer Africa (PTY) Ltd
The Estuaries · Century City
Cape Town, 7441
Phone +27 (0) 21 555 - 1415
Fax +27 (0) 21 555 - 0639
info.southafrica@stieglmeyer.com
www.stieglmeyer.co.za

Nederland / Netherlands
Stieglmeyer BV
Peppelenbos 7
6662 WB Elst
Phone +31 (0) 345 581 - 881
info@stieglmeyer.nl
www.stieglmeyer.nl

Suomi / Finland
Stieglmeyer Oy
Ahjokatu 4 A
15800 Lahti
Phone +358 (0) 3 7536 - 320
Fax +358 (0) 3 7536 - 330
stieglmeyer@stieglmeyer.fi
www.stieglmeyer.fi

Polska / Poland
Stieglmeyer Sp. z o.o.
86-212 Stolno
Grubno 63
Phone +48 56 6771 - 400
Fax +48 56 6771 - 410
info.stolno@stieglmeyer.pl
www.stieglmeyer.pl



Finden Sie Ihren Ansprechpartner auf unserer Website:
Find your local partner on our website:
www.stieglmeyer.com/partner

stieglmeyer.com