



Accessories for hospital beds

Enhanced flexibility and safety



Accessories from Stiegemeyer – Modern solutions for meeting patients' needs

Discover our MDR-compliant customisation options

“Solutions for people who care” is the goal of the Stiegemeyer Group. As a solution provider, we not only look at the hospital bed itself but also at customisation options to meet the needs of patients. With over 170 accessories to choose from, equipping your beds for a wide range of different requirements can be done quickly and tool-free.

And you will benefit from a number of strong advantages. Our accessories comply with the requirements of the EU Medical Device Regulation (MDR) and meet a high standard of safety. Many accessories are suitable for universal use on multiple bed models and thus support flexible and economical working. This is backed up by their long availability – follow-up deliveries and spare parts orders can be supplied as in the customary way. Enhance the care of your patients and relieve the pressure on your staff with the use of original accessories from Stiegemeyer.



**Patient lifting poles
and grab handles**

Page 4



Infusion stands

Page 5



**Mattress base
extension elements**

Page 6



**Accessories on the
footboard**

Page 7



Standard sliding rails

Page 8



Universal holders

Page 10



Oxygen bottle holders

Page 12



**Extension boards and
traction frames**

Page 13



Handset holders

Page 14



Urine bottle holders

Page 15



Disinfectant holders

Page 15



Protection covers

Page 16



Reading lamps

Page 17



Out-of-Bed flex

Page 18



Bed adaptations

Page 19



Overview

Page 20

Patient lifting poles and grab handles

In order to lie comfortably in bed for a lengthy period of time, everyone has to change their position every so often. Our patient lifting pole supports patients in getting themselves into a comfortable position on their own. Stainless steel versions provide particularly good corrosion protection.

The patient lifting pole is inserted into the adapter sleeves on the right or left side of the hospital bed at the head end. The oval or

round metal tube is combined with an ergonomically formed grab handle. The patient can raise himself/herself with little effort. We also offer optional patient lifting poles with double hooks for infusions. They allow safe, space-saving treatment of the patient.

Our cable clips are used to attach handset and push button cables securely to the patient lifting pole.



▲ The oval-shaped stainless steel patient lifting pole, optionally with double hooks for infusions and a soft-touch grab handle.



▲ A special grab handle for the Seta pro junior youth bed provides standard-compliant protection and comfort for the patient.



▲ The zinc-coated cylindrical lifting pole provides stable and reliable service.



▲ Cable clips.

Infusion stands for reliable medication

Infusion stands must be flexibly positionable and offer great stability at the same time. Models attached directly to the bed will not become a stumbling block for patients and visitors and can be positioned close to the patient.

Our infusion stands can be divided into four groups: The chrome-plated standard infusion stand with simple screw adjustment is suitable for most day-to-day tasks. Its hooks can optionally hold 2 or 5 kg.

If you use machine washable beds, we recommend the use of the stainless steel variant with one-hand adjustment.

With its 45° bend, this infusion stand avoids collisions with medical supply systems when the height of the bed is adjusted.

For even higher requirements – such as for rinsing – a heavy-load infusion stand is available with 2 hooks, each carrying 10 kg.



▲ The stainless steel infusion stand is particularly suitable for machine washable beds.



▲ The heavy-duty infusion stand carries loads of up to 10 kg per hook.



▲ The infusion stand with a 45° bend prevents collisions when adjusting the bed.



▲ Alternatively to an infusion stand, a holder with 4 hooks can be attached to the patient lifting pole.


Mattress base extension elements

People in general are growing increasingly taller and consequently need more space in a hospital bed. To make the best use of the integrated extension which our beds offer, we recommend using our tailor-made accessories. In the absence of a pull-out or fold-out

mattress base extension on your bed, you can use our attachable metal insert to extend the lower leg rest without the use of tools. Extend the mattress in every case using our mattress piece to fill the gap between the mattress and the foot end of the bed.



▲ The lower leg rest can be extended fully tool-free with the use of our attachable metal insert. For extending the mattress, additional mattress pieces in 3 heights can be selected.



▲ Some beds can be configured with an integrated extendible lower leg rest extension. Matching mattress extension pieces are also available for this purpose.

Accessories on the footboard for the daily routine

Even in the age of digitalisation, quick handwritten notes and entries are often necessary in everyday clinical life. Our writing surface provides the best conditions for doing so. It is easy to attach to the footboard and can be folded down to save space when not in use. Paper is held securely in place during writing. The large document holder is used to store documents. It also offers space for magazines and is also hung on the handle bars on the headboard or footboard or on the bedside cabinet.

Being able to address patients by name provides an atmosphere of trust and is a help for doctors on their rounds. Our plastic name tag holder for a name tag is practical and hygienic.




▲ The plastic name tag holder can be positioned flexibly.




▲ The writing surface holds the patient file or the tablet PC in place to allow records to be filled in.



▲ The document holder in a modern Argentum colour goes well with every bed.



▲ A handy towel holder can be attached to the handle bars of the headboard and footboard.



▲ The crutch holders can be attached to the bed or bedside cabinet.



Standard-compliant attachment of medical devices according to DIN EN ISO 19054

The hook-on stainless-steel standard sliding rail is a practical solution for attaching intensive care equipment to universal beds.



The adaptable standard sliding rail provides sufficient clearance between the lowered mattress base and the floor.



Secure attachment option for essential accessories: the hook-on short standard sliding rail.



Standard sliding rail on the bed frame (head end).



Standard sliding rail on the backrest.

Standard sliding rails

Standardised attachment of medical products

Stainless steel standard sliding rails are extremely stable and provide a secure hold for medical equipment and drainage bottles on the hospital bed. They are suitable for use in all hospital units, especially in intensive care units.

Hook-on standard sliding rails of different lengths can be attached on all sides of the bed. The bracket on the adaptable standard sliding rail is bent upwards, so that the sliding rail maintains sufficient clearance from the floor even in the lowest position of the bed.

Standard sliding rails are also available for the modern plastic split safety sides in the Evario family and for the Sicuro tera intensive care bed. This means important accessories can be positioned close to the patient and are not concealed even if the safety side is lowered.



Standard sliding rail for the safety side on the Evario one.



Standard sliding rail for the Protega safety side on the Evario.

Accessories for medical products

The EU Medical Device Regulation (MDR) not only formulates precise requirements for hospital beds, but also for accessories.

- Our accessories are specifically designed for use on our beds and provide the best possible **safety and functionality**. Certain general accessories, such as the writing surface or the document holder, can also be used flexibly on other bed models.
- For a medical product, and for accessories for a medical product, a **clinical evaluation** must be carried out.
- Stieglmeyer fulfils its obligation to demonstrate compliance with essential safety and performance requirements as part of a **conformity assessment procedure**.

Our accessories meet all **regulatory requirements** for use on our beds. Therefore, use only original accessories from Stieglmeyer.



The universal accessory rail with plenty of room for accessories.




▲ The mounting bar is suitable for intensive care.



▲ The drainage tube holder can be attached to the handle bars on the headboard and footboard.



▲ Universal holder on a Puro hospital bed with flexible attachment heights for the low-height adjustable mattress base.



▲ The accessory hook can be hung, for example, on various plastic safety sides.

Universal holders for flexible patient care

Accessories, drainage bottles and urine bottles need to be positioned safely along the sides of the hospital bed. Tubes must not be sharply bent or become detached, and liquids must not leak. Use our individual solutions for secure attachment.

External devices such as blood pressure monitors or infusion pumps can be attached to the mounting bar. Various holders can accommodate urine and drainage bags.

To secure additional drainage tubes, practical tube holders can be hung on the headboard or footboard or on the mobilisation support of the MultiFlex safety side.

Versatile attachment options

- Our flexibly attachable holders and accessory bars promote the care of the patient from all sides of the bed.
- Their tool-free installation enables care staff to prepare the bed independently with the minimum of effort.
- Our universal holders are designed specifically for use on Stiegemeyer beds.



The universal accessory holder allows drainage bags and Redon bottles to be attached securely.



On the Puro, oxygen bottles of up to 10 litres can be attached to the holder on the stainless steel standard sliding rail.

Oxygen bottle holders close to the patient

The reliable supply of oxygen to respiratory patients is always a top priority. Oxygen bottles must remain in the required position on the bed, even when the bed is moved. Our oxygen bottle holder enables bottles containing up to 2 litres to be flexibly positioned on the headboard and footboard. The holder is simply hung from

the handle bar and secured with a strap with a Velcro fastener. Its holder fits on all headboard and footboard variants. And the oxygen bottle holder for attaching to the stainless steel standard sliding rail has an even higher safe working load: it holds up to 10-litre bottles.



◀ The stainless steel handle bar on the Rondo footboard provides a secure hold for the oxygen bottle holder.



The foot-end extension fitting on the Evario hospital bed provides ideal conditions for health-promoting therapy of the patient.

Extension boards and traction frames Precise fitting of devices



If extremities need to be fixed in a certain position in trauma surgery units, using perforated rod systems, for example, one of the basic requirements is the ability to attach extension devices safely and precisely to the hospital bed. An exchangeable extension board, which can be attached to the bed in place of the footboard, is used for this purpose. Its two sturdy poles provide an ideal hold for the devices.

◀ A further treatment alternative is the traction frame, which spans the entire bed.

Handset holders

always within easy reach



Handsets should always be within easy reach of patients and care staff and, at the same time, they should be kept in a safe place on the bed to avoid any damage. Our hospital beds already offer practical fixing options as standard.


Handsets with a hook can be hung, for example, from the handle bars of the headboard and footboard. An additional holder can be attached there on request to prevent the handset from slipping or falling off.

We also recommend our handset holder for the purpose of flexibly positioning handsets without hooks. It is ideal, for example, for the innovative LCD handset.

◀ With a handset holder on the headboard or footboard, the LCD handset can be placed in or out of immediate reach of the patient.



▲ There is also a tailor-made holder for the standard handset.



▲ The holder for attaching the handset prevents it from slipping or falling off.

Urine bottle holders

for ideal positioning

Urine bottles should be easy to reach whenever necessary and reliably prevented from tipping over. Our urine bottle holders are ideal for meeting these requirements. We offer models with a holder to suit a variety of hospital beds and equipment options.



▲ The urine bottle holder for the Puro's MultiFlex' safety side with flexible attachment heights for the low-height adjustable mattress base.

The urine bottle holder for attaching to the side of the bed. ▶



▲ The urine bottle holder for the Protega safety side on the Evario.



Disinfectant holders

at the point of care



Reliable hand hygiene is a basic requirement in the fight against nosocomial infections. We offer a holder which allows disinfectant to be positioned right on the bed. This means care staff can disinfect their hands directly before any contact with the patient.

◀ The disinfectant holder promotes reliable hygiene in hospitals.

Protection covers for special requirements

Our protection covers offer enhanced safety for restless patients. They encase full-length and split safety sides and provide protection from impact injuries. They are soft and pleasing to the touch which, in turn, has a soothing effect. In most cases, care staff can raise and lower the safety sides together with the protection

covers. Integrated control panels on the outsides remain accessible. The protection covers are easy to remove for cleaning. Their material is particularly wear-resistant and easy to care for. It can withstand urine, blood and disinfectants.



▲ The protection cover for the full-length safety side moves up and down with it when the bars are raised and lowered.



▲ The split MultiFlex* safety side can also be fully encased in an upholstered cover.



▲ The protection cover pleasantly cushions the four Protega safety side elements.



▲ Lateral tilting of the bed is a recognised means of improving the oxygen supply to intensive care patients. Suitable cushions provide protection and comfort.

Reading lamps for greater comfort for patients

Our convenient reading lamps for attaching to the bed enhance every private room. Every model uses LED bulbs, which offer the advantage of a long life and efficient energy use.

The Sola reading lamp is made of shatterproof plastic. Vents and heat dissipating cylinders ensure that heat build-up remains low.

The Stella reading lamp stands out with its contemporary lamp head made of high-gloss polished stainless steel, as well as with the night light which provides better orientation.

With its cable switch, the attractive Sira reading lamp is particularly easy for the patient to switch on and off. The Sira and Stella models are also available in variants for attaching to the bedside cabinet.



▲ The Sola reading lamp has an attractively shaped lamp head made of shatterproof plastic.



▲ Our Stella reading lamp provides bright lighting – without disturbing the patient in the next bed.



▲ The Sira reading lamp convinces with its compact lamp head and its practical cable switch.


Out-of-Bed flex – the sensor bar for flexible use

With the Out-of-Bed flex mobile assistance system (OOB flex), hospitals are able to cater to the needs of fall-prone patients on a case-by-case basis. The mobile OOB flex system notifies care staff whenever a patient gets up and does not return to bed within a specified period of time. It offers protection and quality of life combined with cost-efficient flexibility. The system can be installed, only when required, on different Stiegmeyer hospital beds. The OOB flex system enables hospitals to meet a tailored investment need, and provide a high standard of care for their patients while also easing the workload of care staff – particularly in situations with fewer people on duty.

The integrated Out-of-Bed system

Alternatively to the OOB flex, our hospital beds can be configured with an integrated Out-of-Bed system right from the outset when you place your order.

- Precise sensor in the backrest prevents false alarms
- Convenient setup and use with an intuitive handset



1 Sensor bar
The delivery includes assembly parts for hospital beds. Usable irrespective of the bed model variant or the mattress base width.



2 Control unit
The control unit can be attached at an ergonomic place to suit the particular care situation on the headboard, footboard or side panels. The power supply is via a separate power plug.

Operating concept
The user interface offers a choice of four different time settings at the touch of a button before a bed exit alert is triggered.




Better orientation at night
A reading lamp can be connected to the automatic light system provided by OOB flex, which switches the light on whenever the patient gets up. This provides better orientation at night and reduces the risk of falls.

Bed adaptations Special solutions for your requirements

Our extensive selection of more than 170 accessories offers a range of combinations to cater for many applications. But special cases requiring a specific solution can still arise – such as in intensive care, in an intermediate care unit or in psychiatric care. In this case, we are pleased to examine the possibility to adapt accessories to your particular requirements. We also offer retrofit

solutions through our service team after assessing and costing your needs.

Please feel free to contact your familiar service or sales team if you would like to enquire about a solution. Together we will find a solution – after all, our slogan is “Solutions for people who care”.



▲ Our extensive portfolio of accessories is suitable for many requirements – in this case, you can see an example of an Evario for intensive care. In special cases, we work together with you to find an individual solution to meet your needs.

Overview

Page		Acute care / Intensive care					
		Evario one	Evario	Sicuro tera	Puro	Seta pro	Deka
4	Patient lifting pole and grab handle	✓	✓	✓	✓	✓	✓
5	Infusion stand	✓	✓	✓	✓	✓	✓
6	Mattress base extension element	✓	✓	✓	✓	✓	✓
7	Accessories on the footboard						
	Name tag holder	✓	✓	✓	✓	✓	✓
	Writing surface	✓	✓	✓	✓	✓	✓
	Document holder	✓	✓	✓	✓	✓	✓
	Towel holder	✓	✓	✓	✓	✓	✓
	Crutch holder	✓	✓	✓	✓	✓	✓
8-9	Standard sliding rail	✓	✓	✓	✓	✓	-
10-11	Universal holder	✓	✓	✓	✓	✓	✓
	Tube holder	-	✓	✓	✓	✓	✓
12	Oxygen bottle holder	✓	✓	✓	✓	✓	✓
13	Extension board	-	✓	✓	✓	✓	✓
14	Handset holder	-	✓	✓	✓	✓	✓
15	Urine bottle holder	✓	✓	✓	✓	✓	✓
15	Disinfectant holder	-	✓	✓	✓	✓	✓
16	Protection cover	✓	✓	✓	✓	✓	✓
17	Reading lamp	✓	✓	✓	✓	✓	✓
18	Out-of-Bed flex	✓	✓	✓	✓	✓	✓

Pediatrics			Transport	
Seta pro junior	Cuddle	Jovie	Mobilo	Ravello
✓	-	-	-	-
✓	✓	-	✓	✓
✓	-	-	-	-
✓	-	✓	-	-
✓	-	-	✓	-
✓	-	-	-	-
✓	-	-	-	-
-	-	-	-	-
✓	✓	-	-	✓
✓	-	-	-	-
✓	-	-	✓	-
-	-	-	-	-
✓	-	-	-	-
✓	-	-	-	-
✓	-	-	✓	-
✓	-	-	-	-
-	-	-	-	-

Please note that not all accessories are compatible across all models. There may be exclusions for certain equipment variants. Please ask your contact person from our sales team for the exact material number which is suitable for your particular product configuration.

Bear in mind that this list of accessories does not claim to be exhaustive. Ask your contact person from our sales team if the accessory you are looking for is not included.

If you have an older product which is not included in the list of compatible models, please ask us about available accessories.

Useful extras

This may also interest you!



▲ With the short or long hook, towels or toiletry bags can be stored close to the patient.



▲ A plug-in USB extension cable allows mobile devices to be charged within easy reach of the patient.



▲ The cup holder provides good support for glasses and bottles.

As a provider of holistic solutions, Stieglmeyer also supplies accessories for other products – such as for bedside cabinets.

There is also a wide variety in the area of mattresses. These are assisted, for example, by modern alternating pressure systems, which provide tailored support for patients on a case-by-case basis.



▲ The Novacare APM 440 XC alternating pressure system is an innovative air cell mattress for immobile patients that offers optimal prevention of pressure ulcers up to stage IV.

STIEGELMEYERGROUP



125 YEARS

EXPERIENCED HUMAN SOLUTION-FOCUSED

STIEGELMEYER GROUP

Solutions for people who care

The Stieglmeyer Group is Germany's leading manufacturer of medical beds for hospitals, nursing homes and homecare. With modern bedside cabinets, furniture and digital assistance systems, it is also an innovative provider of comprehensive solutions. The family-run company successfully develops offers for international markets in more than 100 countries.

Our strengths:

- Product development at our Herford headquarters in Germany – production in four modern factories in Germany and Poland
- Factory in China for the Chinese market
- Subsidiaries and sales partners in many countries
- Top-class service with reliable spare parts supply
- Training and process optimisation for our customers
- Sustainable action with strong social commitment

Our expertise

- Hospital
- Nursing home
- Homecare
- Lifestyle beds
- Furniture systems
- Service

Our values

- Experienced
- Human
- Reliable
- Solution-focused
- Open-minded
- Confident

1 000 EMPLOYEES WORLDWIDE

11 LOCATIONS in Germany, Poland, the Netherlands, Belgium, France, Finland, South Africa and China

Tailor-made OFFERS for all MARKETS

MORE THAN **100 000** beds produced every year

UP TO **15-YEAR** spare part availability

DISTRIBUTION in more than **100** countries

International CERTIFICATION

GERMAN ENGINEERING

Our brands





STIEGELMEYERGROUP

Solutions for people who care



Hauptsitz / Headquarters

Deutschland / Germany
Stieglmeyer GmbH & Co. KG
Ackerstraße 42
32051 Herford
Phone +49 (0) 5221 185-0
Fax +49 (0) 5221 185-100
info@stieglmeyer.com
www.stieglmeyer.com



Hauptsitz / Headquarters

Deutschland / Germany
Burmeier GmbH & Co. KG
Industriestraße 53
32120 Hiddenhausen
Phone +49 (0) 5223 9769 - 0
Fax +49 (0) 5223 9769 - 090
info@burmeier.com
www.burmeier.com

France
Stieglmeyer SAS
16 Avenue de l'Europe
Bâtiment SXB1
67300 Schiltigheim
Phone +33 (0) 4 7238 - 5336
contact@stieglmeyer.fr
www.stieglmeyer.fr

South Africa
Stieglmeyer Africa (PTY) Ltd.
The Estuaries - Century City
Cape Town, 7441
Phone +27 (0) 21 555 - 1415
Fax +27 (0) 21 555 - 0639
info.southafrica@stieglmeyer.com
www.stieglmeyer.co.za

Nederland / Netherlands
Stieglmeyer BV
Peppelenbos 7
6662 WB Elst
Phone +31 (0) 345 581 - 881
info@stieglmeyer.nl
www.stieglmeyer.nl

Suomi / Finland
Stieglmeyer Oy
Ahjokatu 4 A
15800 Lahti
Phone +358 (0) 3 7536 - 320
stieglmeyer@stieglmeyer.fi
www.stieglmeyer.fi

Polska / Poland
Stieglmeyer Sp. z o.o.
86-212 Stolno
Grubno 63
Phone +48 56 6771 - 400
Fax +48 56 6771 - 410
info.stolno@stieglmeyer.pl
www.stieglmeyer.pl

België / Belgique / Belgium
Distrac Group
3320 Hoegaarden
Bleyveld 14
Phone +32 (0) 16 76 82 82
info@distracgroup.com
www.distracgroup.com

China
Stieglmeyer China Co. Ltd.
Poly Center, No. 8 South Wenhua Road
Dongping New Town, Shunde District
528300 Foshan
Phone +86 757 29878710 - 0
info@stieglmeyer.cn
www.stieglmeyer.cn



stieglmeyer-group.com