

Bedside cabinets and tables



Vitano

Flexible use with rotatable body

- Cabinet body rotates in both directions
- With additional electrical options
- Overbed table has a high load-bearing capacity
- Closed back prevents unauthorised access from behind



Somero

The bedside cabinet with or without an overbed table

- Easy-to-clean fronts
- Overbed table easy to change from one side to the other
- Closed back prevents unauthorised access from behind



Soreno eco

The bedside cabinet with space-saving storage

- Two drawers in different sizes
- Closed back prevents unauthorised access from behind
- Budget-oriented model offers excellent value for money



Quado

More than a bedside cabinet

- Height adjustable from 85 to 115 cm
- The overbed table folds down easily and softly
- Drawers open from both sides
- Machine washable
- Also ideal as a multifunctional trolley, ward trolley or cleaning trolley

Mattresses



▲ Z1: Dimensions (WxLxH): 87/90 x 200 x 15 cm



▲ Z2: Dimensions (WxLxH): 87/90 x 200 x 12 cm



▲ Dimensions (WxLxH): 87/90 x 200 x 15 cm



▲ P1: Dimensions (WxLxH): 87/90 x 200 x 16 cm

P2: Dimensions (WxLxH): 87/90 x 200 x 14 cm ▲



▲ C1: Dimensions (WxLxH): 87/90 x 200 x 12 cm



▲ S1: Dimensions (WxLxH): 87/90 x 200 x 16 cm



Z1 – Sky-therm®

with wavy cut

- Wavy cut protects problem areas of the body
- Excellent point elasticity due to lattice of the Sky-therm® layer

Z2 – Soft-therm

with fine perforations

- Highly elastic soft-therm layer
- Fine perforations enhance the flexibility and improve the microclimate

TEMPUR-MED®

PRESSURE REDISTRIBUTING MATTRESSES AND PILLOWS

Premium mattress

- Outstanding pressure ulcer prevention validated in trials
- Reacts to body weight and temperature

P1, P2 – Prevention mattresses

Protection against decubitus

- Strong visco-foam layer ensures relieving pressure distribution
- Edge zone reinforcement for greater protection and mobility

C1 – Comfort mattress

Great reclining properties

- „Sandwich cushion“ made of visco and cold foam with good microclimate

S1 – Heavyweight mattress

for high requirements

- 3-layer mattress for heavier users

Alternating pressure systems for clinical use

- Systems for different requirements
- Durable PU membranes
- Acoustic signal functions
- Control units with robust housings and retaining brackets



Subject to technical changes without notice. For technical reasons, the actual colours may differ slightly from those shown in print. Last updated 01/2024, Version 12, Item no. 601574.

STIEGELMEYERGROUP



Our products for hospitals

Beds, bedside cabinets and accessories

STIEGELMEYERGROUP
Solutions for people who care



Hauptsitz / Headquarters

Deutschland / Germany
Stiegelmeier GmbH & Co. KG
Ackerstraße 42
32051 Herford
Phone +49 (0) 5221 185-0
Fax +49 (0) 5221 185-100
info@stiegelmeier.com
www.stiegelmeier.com

stiegelmeier-group.com



Hauptsitz / Headquarters

Deutschland / Germany
Burmeier GmbH & Co. KG
Industriestraße 53
32120 Hiddenhausen
Phone +49 (0) 5223 9769-0
Fax +49 (0) 5223 9769-090
info@burmeier.com
www.burmeier.com

International Edition

ESSENTIALLY BETTER



High quality medical products from Stiegmeyer

Beds and furniture for use in hospitals must cover a great spectrum of requirements. Safety in everyday use is the key consideration. Outstanding hygiene properties are important to win the battle against drug-resistant superbug infections. Patients' recovery should be supported through high comfort. Nursing staff rely on easy ergonomic bed controls.

Stiegmeyer has been successfully working to meet these requirements for many decades and to provide value for money products. The family firm from Herford has long been a pioneer of many core areas of bed development, and still is today. The design department can look back on more than forty years' experience with machine washable hospital beds.

The Stiegmeyer catalogue contains beds for both special mobilisation or intensive care requirements as well as models offering an entire modular system for customisable configurations. This makes it possible to meet the individual needs of every institution down to the last detail. The quality of Stiegmeyer is characterised by lifelong support and 15 years spare part availability. Further service offers can be requested from our local partners.

Hospital beds from Stiegmeyer introduce a new feeling of well-being into the daily routine of hospital units.

Hospital beds



Evario one

with split 50/50 safety sides

- Height adjustment range: approx. 36–82 cm
- Safe working load 250 kg
- 50/50 safety side offers full protection
- Optional with fixed headboard and radiolucent backrest
- Operation using the handset, integrated control panel or nurse control



Evario one

with split 3/4 side safety side

- Height adjustment range: approx. 36–82 cm
- Safe working load 250 kg
- The head end of the safety side can support mobilisation
- Operation using the handset, integrated control panel or nurse control



Evario one

with pivoting 3/4 safety sides

- Height adjustment range: approx. 36–82 cm
- Safe working load 250 kg
- Pivoting safety sides
- Operation using the handset or nurse control
- Attached accessories are always visible



Evario

with Protega safety sides

- Height adjustment range: approx. 35–91 cm
- Safe working load 260 kg
- Safety sides with integrated bed control
- Optional integrated bed scale with an accuracy of 200 g



Evario

with 3/4 side safety side

- Height adjustment range: approx. 32–91 cm
- Safe working load 260 kg
- Innovative LCD handset with only 3 buttons
- Quiet lowering of the safety side

Hospital beds, transport, stretchers



Sicuro tera

The bed for intensive care

- Height adjustment range: approx. 45–84 cm
- Safe working load 260 kg
- Lateral tilting to both sides up to 25°
- Optional integrated bed scale with an accuracy of 200 g



Puro

Safe and hygienic

- Height adjustment range: approx. 28–91 cm
- Safe working load 275 kg
- MultiFlex+ safety side
- Excellent machine washability
- LCD handset with only 3 keys and a backlit display



Seta pro

The customisable bed

- Height adjustment range: approx. 32–91 cm
- Safe working load 260 kg
- Full-length safety side integrated in head/footboards
- Optional narrow bed with an external width of less than 100 cm



Deka

The modular platform

- Height adjustment range: approx. 43–83 cm
- Safe working load 225 kg
- Mechanical and motor versions
- Optionally machine-washable
- Also available as a psychiatric bed



Mobilo

The stretcher for transporting

- Safe working load 200 kg
- 2-section mattress base (optionally also 4-section)
- Safety sides optionally on both sides
- High comfort

Children's beds and accessories



Jovie

The modern baby cot

- Lying surface height: approx. 50 - 75 cm
- Suitable for babies with a height up to 60 cm and a weight of up to 12 kg
- Safe working load 20 kg
- Foot rest position up to 13°
- Compliant with EN 50637



Cuddle

The friendly children's bed

- Height adjustment range: approx. 77–111 cm
- Safe working load 100 kg
- Mattress base 72.5 x 132.6 cm
- Suitable for children with a body height up to 125 cm
- Compliant with EN 50637



Seta pro junior

The hospital bed for young patients

- Height adjustment range: approx. 35–91 cm
- Mattress base 78 x 180 cm
- Suitable for children with a body height up to 155 cm
- Full-length safety side integrated in head/footboards
- Compliant with EN 50637



Ravello

The transport chair

- Safe working load 200 kg (static), 135 kg (dynamic)
- 43 degree tilting backrest
- Lowerable arm rests
- Headrest



Accessories

for hospital beds

- Patient lifting poles and infusion stands
- Urine bottle and oxygen baskets
- Refer to our accessory catalogue for more information